



Stereo Microscope Setup & Use

1. Introduction

A stereo microscope is an essential tool for examining insects in detail. It provides:

- A **three-dimensional view** of specimens
- Adjustable magnification for **general observation and detailed study**
- Safe examination of small or fragile specimens without damaging them

Common uses in entomology:

- Sorting and identifying field-collected specimens
- Observing structural features (legs, antennae, wings, setae)
- Preparing and checking pinned specimens
- Dissecting small body parts (e.g., genitalia, wings)

2. Setting Up the Microscope

1. Workspace

- Clean, flat, stable surface free from clutter
- Good lighting in the room; avoid strong backlight reflections
- Keep water, drinks, and food away

2. Positioning the Microscope

- Place the microscope at **comfortable eye level**
- Adjust height to avoid leaning forward excessively

3. Power and Illumination

- Plug in or ensure batteries are charged (if LED)
- Choose **reflected light** for opaque specimens (beetles, wings)

- Use **transmitted light** for translucent parts (wings, thin appendages)
- Adjust intensity gradually; avoid overly bright light

4. Eyepiece Adjustment

- Adjust **interpupillary distance** so both eyes see a single image
- If your microscope has **diopter adjustment**, set it according to your dominant eye

3. Magnification Basics

- Stereo microscopes typically have **fixed zoom (e.g., 10×–40×)** or **rotating magnification turret**
- Start at **lowest magnification** to locate your specimen
- Gradually increase magnification to examine fine details
- Use **coarse focus** first, then fine focus for clarity

Tip: Keep notes on **magnification used** for reference in your field or lab records.

4. Handling Specimens Safely

- Use **soft brushes or fine forceps** to move specimens
- Place small insects on **white or black trays** for contrast
- Pin specimens **carefully**; avoid forcing them under the lens
- Keep specimens **stable** – a piece of foam or modelling clay under a tray helps

5. Observation Techniques

- **Avoid eye strain:**
 - Take breaks every 15–20 minutes
 - Adjust chair and microscope height for posture
- **Adjust lighting:**
 - Angle LED lights to highlight textures and shadows
 - Use both reflected and transmitted light for different features
- Use **low magnification to orient**, then zoom in for detail

- **Record observations:**
 - Sketch or photograph structures
 - Note magnification, specimen ID, and date

6. Practical Hacks & Tips

1. **Use a ruler under the specimen** for quick size estimation
2. **Pin duplicate specimens** for comparison without moving originals
3. **Small pieces of black paper** can improve contrast for light-coloured insects
4. **Turn off lights between specimens** to reduce heat exposure
5. **Keep a magnifying hand lens nearby** – sometimes faster for rough sorting
6. **Label specimens immediately** to avoid confusion

7. Cleaning & Care

- Clean **eyepieces and objective lenses** with lens tissue only
- Keep **dust covers on** when not in use
- Avoid touching lenses with fingers
- Store microscopes on a **stable surface**; don't move while in use

8. Summary

Using a stereo microscope properly allows:

- Careful observation of insect morphology
- Safe handling and preparation of specimens
- Accurate field and lab records

Remember: Comfortable posture, proper lighting, and gentle handling of specimens will improve both your observations and the longevity of your equipment.

